



Sotair®

- Sotair is a flow control valve that fits between a manual resuscitator and a patient mask/airway. The valve prevents flow rates of air from exceeding 55 LPM during manual ventilation
- Regulation of airflow helps reduce peak airway pressures, excessive tidal volumes¹ and excessive breath delivery rates,² minimizing the risk of over-pressurization, over-ventilation and hyperventilation
- Sotair acts as a 'Forcing Function,' constraining the provider to deliver safer manual ventilation with less variability. This is most effective in ensuring system-wide change^{1 3}
- Immediate, real time haptic, auditory and visual feedback simplifies training
- 100% activation and flow limiting testing ensures reliability in performance
- Adult, single patient use



Aggregate Improvement in Performance



Of Providers Recommend Sotair



Sotair Training Time




1. Brady MF, Weber NK, Walker, III R, et al. Feasibility of manual ventilation replacing mechanical ventilation. BMJ Innovations 2021;7:297-301.
 2. Memphis Sotair Limited Market Release, 2022.
 3. <https://www.ncbi.nlm.nih.gov/books/NBK225188/>, <https://ismpcanada.ca/resource/hierarchy-of-effectiveness/>



Learn Sotair[®] in 30 Seconds

1. **Squeeze the bag fast and forcefully to activate Sotair's ALARMS.** You will *feel* the bag immediately become rigid, *hear* the sound of the Sotair's valve collapsing, and *see* the valve close.

2. **Gradually decrease the speed of your squeeze until Sotair's alarms do not activate.** If the alarms keep activating, slow down your squeeze further until no alarms are experienced



3. **Find the sweet spot: Squeeze over ~1 second, near Sotair's alarm activation, but avoid activation of alarms.**

***Squeeze the bag with 1 hand, over 1 second, every 6 seconds**