

**Four things in one pack,
one less thing to think about**

 **INTERSURGICAL**
COMPLETE RESPIRATORY SYSTEMS
i-gel[®] O₂

Resus Pack

Stay focused on what matters

 **INTERSURGICAL**
COMPLETE RESPIRATORY SYSTEMS

Quality, innovation and choice



www.i-gel.com

INTERSURGICAL
COMPLETE RESPIRATORY SYSTEMS

i-gel^{O₂}

Resus Pack



In emergency medicine, you need equipment that's easy, rapid and reliable to use. The i-gel O₂ Resus Pack contains everything you need to prepare, insert and secure the i-gel O₂ quickly and efficiently: an i-gel O₂ supraglottic airway, a sachet of lubricant, airway support strap and a suction tube.

The i-gel O₂ incorporates a supplementary oxygen port and an integral hook ring for securing of the airway support strap.


The i-gel O₂ has been designed to facilitate ventilation as part of standard resuscitation protocols, such as those designated by the European Resuscitation Council (ERC) and the American Heart Association (AHA). However, the i-gel O₂ incorporates a supplementary oxygen port, so can also be used for the delivery of passive oxygenation, or Passive Airway Management (PAM™), as part of an appropriate CardioCerebral Resuscitation (CCR) protocol.



The i-gel O₂ Resus Pack is available in 3 adult sizes and includes:

- i-gel O₂ supraglottic airway
- Sachet of lubricant – for quick and easy lubrication of the i-gel O₂ prior to insertion
- Airway support strap – to secure the i-gel O₂ in position
- 12 FG suction tube – for insertion through the gastric channel to empty fluid from the stomach.

The new i-gel O₂



Hook ring

To secure the i-gel O₂ in position with the airway support strap.

Gastric channel

Enhances patient safety by providing a mechanism for the management of regurgitant fluid.

Non-inflatable cuff

Eliminates the need for cuff inflation after insertion, allowing easy and rapid insertion.

Supplementary oxygen port

For the administration of passive oxygenation as a component of cardiocerebral resuscitation (CCR).

Integral bite block

Reduces the possibility of airway channel occlusion.

Buccal cavity stabiliser

Aids insertion and eliminates the potential for rotation.

Epiglottic rest

Reduces the possibility of epiglottis 'down folding' and airway obstruction.

Airway innovation since 2007

i-gel is the innovative second generation supraglottic airway device from Intersurgical launched in 2007. The first major development since the laryngeal mask airway, the i-gel changed the face of airway management and is now widely used in anaesthesia and resuscitation around the world.



The non-inflatable cuff


The i-gel O₂ gets its name from the innovative soft, gel-like material from which it is made. It is the innovative application of this material that has enabled the development of a unique non-inflatable cuff. This means there is no need for cuff deflation prior to insertion and no cuff inflation after placement to secure a seal, shortening and simplifying the preparation and insertion procedure.

Insertion

The i-gel O₂ is incredibly easy to use. Insertion is rapid and can normally be achieved in less than 5 seconds. The simplicity of i-gel means the insertion technique is instinctive, reliable and does not require or benefit from the insertion of fingers or thumbs into the patient's mouth.

Fixation




The i-gel O₂ Resus Pack includes a specially designed airway support strap for securing the i-gel O₂ in position. Fixation is easily achieved by sliding the strap under the patient's neck and securing it in place to the hook ring on the collar of the i-gel O₂. This makes it ideal for use where adhesive tape is unsuitable, such as on patients with facial hair.


INTERSURGICAL
 COMPLETE RESPIRATORY SYSTEMS
i-gel^{O₂}
Resus Pack

Ordering information

The i-gel O₂ Resus Pack contains an i-gel O₂, with an easily removable inner tray, containing a sachet of lubricant, airway support strap and suction tube, all packaged in a specially designed clear, sterile, rigid pack. This unique design ensures protection of the contents, both in transit and in storage. The i-gel O₂ has a colour coded hook ring to allow quick and easy identification of size. The i-gel O₂ Resus Pack is available in three adult sizes.



| Code | Description | Box Qty |
|---|--|----------|
|  8705000 | i-gel O ₂ Resus Pack, large adult – includes a size 5 i-gel O ₂ with orange hook ring, sachet of lubricant, airway support strap and a 12FG suction tube (90+kg) | 6 |
|  8704000 | i-gel O ₂ Resus Pack, medium adult – includes a size 4 i-gel O ₂ with green hook ring, sachet of lubricant, airway support strap and a 12FG suction tube (50-90kg) | 6 |
|  8703000 | i-gel O ₂ Resus Pack, small adult – includes a size 3 i-gel O ₂ with yellow hook ring, sachet of lubricant, airway support strap and a 12FG suction tube (30-60kg) | 6 |



Interact with us



Visit the i-gel web site www.i-gel.com

IS6.6 - DC Issue 2 11.13



Intersurgical Ltd,
 Crane House, Molly Millars Lane,
 Wokingham, Berkshire, RG41 2RZ, UK
 T: +44 (0)118 9656 300 F: +44 (0)118 9656 356
 info@intersurgical.com www.intersurgical.com

Deutschland
 info@intersurgical.de

France
 info@intersurgical.fr

Italia
 info@intersurgical.it

España
 info@intersurgical-es.com

Portugal
 info@intersurgical.pt

Nederland
 info@intersurgical.nl

Lietuva
 info@intersurgical.lt

Россия
 info@intersurgical.ru

Česká Republika
 info@intersurgical.cz

Taiwan
 info@intersurgical.com.tw

South Africa
 info@intersurgical.co.za

Philippines
 info@intersurgical.ph

Japan
 info@intersurgical.co.jp

USA
 support@intersurgicalinc.com

China
 info@intersurgical-cn.com

