

AED FEATURE COMPARISON CHART

RESCUE ESSENTIALS

info@rescue-essentials.com
3811 International Blvd., Ste 100, Leland, NC 28451
910-830-0286



	CARDIAC SCIENCE	DEFIBTECH		HEARTSINE		PHYSIO CONTROL	ZOLL	
Model Name	Powerheart G5	Lifeline	Lifeline View	Samaritan PAD 350P/360P	Samaritan PAD 450P	Lifepak 1000	AED PLUS	AED 3
RE SKU	56-0201 (Semi-Automatic) 56-0200 (Fully Automatic)	56-0002	56-0001	56-0210 (350P) 56-0211 (360P) 56-0212 (Connected 360P)	56-0220	56-0000	56-0100	56-0101 (Semi-Automatic) 56-0102 (Fully-Automatic)
Primary Customer	EMS/ Fire	LE/ Public Access	LE/ Public Access	LE/ Public Access	LE/ Public Access	Fire/ EMS	Fire/EMS/ Public Access	Fire/ EMS/ Public Access
Warranty	8 years	8 years	8 years	8 years	8 years	5 Years	7 years	8 years
Weight (lbs)	5.7	4.2	<3	2.4	2.4	7.1	6.7	5.5
Dimensions (H x W x D)	11.8" x 9" x 3.4"	11.8" x 8.5" x 2.7"	9.5" x 7.3" x 2.3"	8.0" x 7.25" x 1.9"	8.0" x 7.25" x 1.9"	3.4" x 9.2" W x 10.9" D	11.5" x 9.5" x 5.25"	5" x 9.7" x 9.3"
Fully-Automatic Option	☑	☑	--	☑	--	☑	☑	☑
Joule Output	95J - 354J	150J	150J	150J - 200J	150J - 200J	150-360J	150J - 200J	150J - 200J
Ingress Protection (resistance to dust, water, vibrations, etc.)	55	54	55	56	56	55	55	55
Data Download Method	USB	Data Cord	Data Cord or USB	USB	USB	USB	USB	WiFi or USB
Electrode Pads Life	2 years	2 years	2 years	4 years	4 years	4 years	5 years	5 years
Electrode Pads Self Test	☑	--	☑	☑	☑	--	--	--
Pediatric Method	Pediatric Pad	Pediatric Pad	Pediatric Pad	Pediatric Pad	Pediatric Pad	Pediatric Pad	Pediatric Pad	Button
Battery Standby Life	4 years	5 or 7 years	4 years	4 years	4 years	4 Years	5 years	5 years
CPR Coaching	Responsive	Passive	Passive	Passive	Responsive	Responsive	Responsive	Responsive
Bilingual	☑	--	--	--	--	☑	--	--
Remote Monitoring	--	--	--	--	--	--	--	☑
Aviation Option	--	--	☑	☑	--	--	--	--
ECG Display	--	--	Optional	--	--	☑	Optional	☑
Convert AED to Trainer	--	☑	☑	--	--	--	--	--
Stand-Alone Trainer (Sold Separately)	☑	☑	--	☑	☑	☑	☑	☑