

The **Fastest, Safest, Most Effective** Prehospital Field Tourniquet

GA7

the NEW **Gen 7**

CAT® COMBAT APPLICATION
TOURNIQUET®



1 **Single Routing Buckle**



- Faster application
- Effective slack removal
- Decreased blood loss
- Fewer windlass turns
- Simplified training with single protocol directions for all applications



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

LABORATORY TESTED COMBAT PROVEN®

PRODUCTS WITH A MISSION®

The **Fastest, Safest, Most Effective** Prehospital Field Tourniquet

GA7

the NEW **Gen 7**

CAT® COMBAT APPLICATION
TOURNIQUET®



2 Windlass Rod



- Increased diameter
- Enhanced strength
- Aggressive ribbing
- Improved grip



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

LABORATORY TESTED COMBAT PROVEN®

PRODUCTS WITH A MISSION®

The **Fastest, Safest, Most Effective** Prehospital Field Tourniquet

GA7

the NEW **Gen 7**

CAT® COMBAT APPLICATION
TOURNIQUET®



3 Windlass Clip

- Rapid windlass lock
- Bilateral buttress
- Added strength



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

LABORATORY TESTED COMBAT PROVEN®

PRODUCTS WITH A MISSION®

The **Fastest, Safest, Most Effective** Prehospital Field Tourniquet

GA7

the NEW **Gen 7**

CAT® COMBAT APPLICATION
TOURNIQUET®



4 Windlass Strap



- Sonic welded to clip for constant contact
- Color changed to Gray for tactical considerations (Black C-A-T® only)



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

LABORATORY TESTED COMBAT PROVEN®

PRODUCTS WITH A MISSION®

The **Fastest, Safest, Most Effective** Prehospital Field Tourniquet

GA7

the NEW **Gen 7**

CAT® COMBAT APPLICATION
TOURNIQUET®



5

5 Stabilization Plate

- Reinforced, beveled contact bar



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

LABORATORY TESTED COMBAT PROVEN®

PRODUCTS WITH A MISSION®

The **Fastest, Safest, Most Effective** Prehospital Field Tourniquet

GA7

the NEW **Gen 7**

CAT® COMBAT APPLICATION
TOURNIQUET®



6 Red Tip Technology

- Red elliptical tip
- Added visual cue for application

7 Free-Moving Internal Band

- Patented band within a band
- True evenly distributed circumferential pressure



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

LABORATORY TESTED COMBAT PROVEN®

PRODUCTS WITH A MISSION®