

DELTA GATEWAY

ENABLING APRU 6L CELLULAR CONNECTION...ANYWHERE



BENEFITS

Connects your APRU 6L portable blood refrigerator to the Cloud (via Bluetooth)

Access temperature data through BloodCOMM app or web browser

APRU 6L cold chain data accessible from anywhere in the world (BloodCOMM subscription required)

Transmits temperature logs to the cloud for constant monitoring of blood or temperature sensitive medications and vaccines

Enables remote alerts when the APRU 6L battery is low or approaching a temperature excursion

Mounting: tabletop, wall mount, or Velcro onto the APRU 6L

AT&T and Verizon Carriers

SIM card included and automatically activated with purchase

SPECIFICATIONS

Part Number: A-001732

Dimensions:

6" x 4" x 5.5" (antenna extended vert.)

10.5" x 4" x 1.5" (antenna extended horiz.)

Weight: 170 grams (0.37 lbs)

Maximum distance of 15 feet from APRU

Temperature rating: -40° to +80° C (on DC power)

Water resistant:

Basic package for indoor use only

IP65 upgrade available

