

RE SIM Wound Packing Task Trainer



Affordable Wound Packing Training

RE SIM from rescue Essentials brings practical, budget-friendly realism to a wide range of instruction, from entry-level civilian courses to advanced Tactical Emergency Casualty Care (TECC) and Tactical Combat Casualty Care (TCCC) programs used by first responders, combat medics, and corpsmen.

The RE SIM Wound Packing Task Trainers are designed to support high-frequency training in wound packing and the application of hemostatic gauze without driving up costs. Each model replicates realistic external anatomy and includes wound cavities with hidden bleed sources to simulate active hemorrhage. Trainers are available in two configurations: a laceration or puncture-style wound, and a large-caliber gunshot wound (GSW). The GSW model incorporates a simulated bone structure to reinforce proper technique when training for full-depth packing to the bone.

Built from rugged silicone, these trainers are made for repeated use and do not require batteries, external power, or mechanical components. A reinforced connection point helps keep tubing secure during training, while an integrated one-way valve prevents backflow of simulated blood into the squeeze bottle. After use, the trainer can be rinsed clean with water, even if the simulant has dried.

Fundamentals of Wound Packing

- Clearing and assessing the wound
- Identify the source of the bleeding
- Preparing wound packing while hands are slick or bloody
- Packing deep to the bone, targeting the source of bleeding
- Keeping firm, consistent pressure

RE SIM Task Trainers offer a reliable way to prepare students for real-world hemorrhage control across TCCC, TECC, TEMS, PHTLS, and other bleeding control and advanced gauze application scenarios.

SPECIFICATIONS

Gunshot Wound:

REF 30-0920V2-D

REF 30-0920V2-L

Laceration:

REF 30-0928V2-D

REF 30-0928V2-L

Dimensions:

2.75" H x 5" W x 8" L
(7cm H x 13cm W x 20cm L)

Weight:

GSW : 2 lbs 7 oz (1.1 kg);
Laceration: 2 lbs 5 oz (1 kg)

FEATURES

- Available with large-caliber gunshot wound with bone or laceration/stab wound
- Enables continual practice to strengthen muscle memory
- Simulated bone included in GSW version for accurate sensory wound packing
- Comes in dark or light flesh tones
- Made from durable silicone
- May be packed dry or with simulated blood
- Easy cleanup and storage
- Includes an additional one-way valve