

PRICING

2022 TRUECLOT® CATALOG



www.rescue-essentials.com

CAPABILITIES STATEMENT

CUSTOM-BUILT TRAUMA AND MEDICAL KITS

Rescue Essentials® is committed to sourcing and delivering custom-built tactical emergency medical kits that improve outcomes in traumatic medical situations. The broad product line and production facility is configured for large runs, and in-house prototyping.

CORE COMPETENCIES

- A complete offering of high-quality, need-driven, tactical situation tested equipment
- Exceptional customer service, with knowledgeable, experienced and dedicated on-site support team
- Specialty prototyping for meeting specific custom specifications
- Production and logistics infrastructure ensuring on-time worldwide delivery
- Single SKU for custom kits simplifies buying process
- Contract assembly for customer specified components
- 10+ years of customer success

BROAD LINE FROM TOP VENDORS

- Sourcing and custom-configuration from wide line of industry leading brands

UNWAVERING COMMITMENT

- A service-first mindset was a Rescue Essentials® hallmark from the beginning, and remains so to this day. We make every effort to stay connected with our customers. This culture of service is shared by all employees and is a key to our success.
- We love what we do and are proud to produce great work, continuously striving to improve outcomes for our customers.

TRUECLOT® PRICING

Item #	Description	MSRP
30-0422	Blood Simulant, Pre-Mix, 1/2 gallon	\$30.00
30-0454	Tourniquet Application Trainer	\$645.00
30-0578	Blood Simulant, Concentrate, 1 gallon + Bottle	\$50.00
30-0662	Wound Packing Training Kit - Lac, Rolled Gauze	\$650.00
30-0670	Inguinal Trainer	\$645.00
30-0771	Wound Packing Training Kit - GSW, Rolled Gauze	\$650.00
30-0777	Basic Packing Trainer (BPT)	\$135.00
30-0912	PCT3 Leg Trainer	\$495.00
30-0916	Blood Simulant, Concentrate, 1 gallon	\$45.00
30-0917	Replacement 1 L Squeeze Bottle/Tubing & Valve	\$29.99
30-0918	Wound Packing Training Kit - GSW, Z-Folded Gauze	\$650.00
30-0919	Blood Simulant, Concentrate, 4 x 1 L Refills	\$60.00
30-0920	Task Trainer, Gunshot Wound w/Bone Insert	\$475.00
30-0921	Simulated Hemostatic Gauze, 12' Z-Folded, 3-Pack	\$25.00
30-0923	Wound Packing Training Kit - Lac, Z-Folded Gauze	\$650.00
30-0926	Bleeding Control Instructors Kit - 1 GSW/1 Knife Wound	\$1,450.00
30-0928	Task Trainer, Laceration/Stab Wound	\$475.00
30-0929	Replacement One-Way Valves for Task Trainers, 5	\$34.00
30-0933	Blood Simulant, Pre-Mix, 1 gallon	\$50.00
30-0936	Bleeding Control Instructors Kit - 2 Gun Shot Wounds	\$1,450.00
30-0942	Bleeding Control Instructors Kit - 2 Knife Wounds	\$1,450.00
30-0943	Bleeding Control Instructors Kit - w/Tourniquet Trainer	\$1,475.00
30-0991	Basic Packing Trainer Instructors Kit	\$1,585.00
40-0318	Simulated Hemostatic Gauze, 4' Z-Folded	\$5.50
40-0321	Simulated Hemostatic Gauze, 12' Rolled, 3-Pack	\$25.00
40-3001	Clotting Solution, 500 mL	\$30.00
61-0067	Replacement 1 L Mixing Bottle	\$15.00
70-0668	Replacement Black Case for Task Trainer	\$35.00



Bleeding Control Instructors Kit with Tourniquet Trainer

The intuitive TrueClot® Task Trainer allows wet or dry wound packing training. The simulated large caliber gunshot wound can be packed with Simulated Hemostatic Gauze while TrueClot® Blood Simulant is pumped into the wound using the attached 1 liter squeeze bottle. The TrueClot® Tourniquet Trainer demonstrates realistic need for proper placement, application and pressure required for efficient hemorrhage control.

Dimensions: 16"H x 20"W x 8"D
Weight: 20 lbs

30-0943 Price: \$1,475.00



Bleeding Control Instructors Kit

Based on the innovative TrueClot® Task Trainers, this kit can be ordered with your choice of GSW or Knife/Laceration wounds. The TrueClot® Task Trainers are designed to simulate wound packing and demonstrate the use of advanced hemostatic gauzes with TrueClot® Blood Simulant. Either wound can be packed with Simulated Hemostatic Gauze while simulant is pumped into the wound using the attached 1 liter squeeze bottle.

Dimensions: 17"H x 20"W x 7.5"D
Weight: 14 lbs

30-0926 1 GSW/1 Knife Price: \$1,450.00
30-0936 2 GSW
30-0942 2 Knife Wounds



Wound Packing Training Kit

This kit includes your choice of either a Simulated Large Caliber GSW with Bone trainer or Knife/Laceration Wound trainer. Plus 4 TrueClot® Blood Simulant Concentrate packs, bottles for mixing and delivering simulant and 9 packs of Simulated Hemostatic Gauze.

Dimensions: 14"H x 10.5"W x 6"D
Weight: 5 lbs

30-0771 GSW/Rolled Gauze Price: \$650.00
30-0918 GSW Z-Folded Gauze
30-0662 Laceration/Rolled Gauze
30-0923 Laceration/Z-Folded Gauze



Basic Packing Trainer Instructors Kit

This kit is based on the innovative TrueClot® Task Trainer (1 GSW trainer) and the unique Basic Packing Trainer (BPT) (6 training kits included). The BPT can be used with or without the included blood simulant. This kit is an excellent for STB, TCCC, TECC, TEMS, and PHTLS courses.

Dimensions: 17"H x 20"W x 7.5"D
Weight: 27 lbs 5.6 oz

30-0991 Price: \$1,585.00



Replacement Components

30-0917
1 L squeeze bottle/tubing and valve
Price: \$29.99

30-1401
1 L mixing bottle
Price: \$12.00

70-0668
Black Case for Task Trainer
Price: \$35.00

30-0929
Replacement One-way Valves for Task Trainers, 5
Price: \$34.00



Tourniquet Trainer

This wearable device allows full application of a tourniquet for training purposes. Training blood can be pumped through the device to simulate hemorrhage, and when a tourniquet is correctly applied, bleeding is stopped. The padded arm cuff reduces pinching and pain and allows for proper application of a tourniquet during classroom or scenario training.

Dimensions: 12"H x 16"W x 7"D (packaged)
Weight: 5 lbs

30-0454 Price: \$645.00



Inguinal Trainer

This wearable training aid simulates a junctional hemorrhage that cannot be treated with a tourniquet and teaches advanced wound packing skills. The trainer can be placed on a manikin or worn by an actor for a higher level of realism. Training blood is pumped through the device to simulate hemorrhage and the large gunshot wound cavity can be packed with a full 12' x 3" strip of simulated hemostatic dressing.

30-0670 Price: \$645.00



Wound Packing Task Trainer

Choose between a GSW or Knife Laceration wound. The Task Trainer features realistic internal and external wound anatomy and ships with 1 complimentary liter of TrueClot® Blood Simulant concentrate and 5 complimentary packs of 4' Simulated Hemostatic training gauze.

Dimensions: 6.5"L x 5.5"W x 2"H
Liquid Capacity: 1 L
Weight: 3 lbs

30-0920 (GSW) Price: \$475.00
30-0928 (Knife/Laceration)



PCT3 Leg Trainer

The TrueClot® Packing, Compression, and Tourniquet Task Trainer (PCT3) is a low cost, durable bleeding control training aid. The simulated leg is ideal for introducing and reinforcing techniques including wound packing, compression bandages, and tourniquet application for hemorrhage control. Tourniquet sold separately.

Dimensions: 13"H x 7" Diameter (at thigh)
Weight: 6.2 lbs

30-0912 Price: \$495.00



Blood Simulant & Solution

Realistic simulated clotting blood for use in TrueClot® Hemorrhage Task Trainers, moulage and other trauma simulation scenarios.

30-0422 (Half-gallon pre-mixed) Price: \$30.00

30-0933 (1 gallon pre-mixed) Price: \$50.00

30-0919 (concentrate, 4x1 L refill) Price: \$60.00

30-0916 (concentrate, 1 gallon) Price: \$45.00

30-0578 (concentrate, 1 gallon + bottle) Price: \$50.00

40-3001 (clotting solution, 500 mL) Price: \$30.00



Basic Packing Trainer

The BPT is a low-cost, introductory level dry wound packing task trainer for the "first exposure" stage of training. Students can repeatedly practice correct wound packing techniques and develop muscle memory to hone their bleeding control skills. The simulated GSW can be packed with 12' of gauze with heavy external pressure. Ideal for table top classroom use.

Dimensions: 2.5"H x 4" Diameter (packaged)
Weight: 13 oz

30-0777 Price: \$135.00



Simulated Hemostatic Gauze

TrueClot® Simulated Hemostatic Gauze works in conjunction with our TrueClot® Blood Simulant to create realistic simulated blood clots.

40-0321 (12' Rolled, 3-Pack) Price: \$25.00

30-0921 (12' Z-Folded, 3-Pack) Price: \$25.00

40-0318 (4' Z-Folded) Price: \$5.50

