

Always Better®



Lightweight



The Medium Trauma Kit NOW!™ (TKN™) replaces a typical Individual First Aid Kit with its unique ability to deploy with one hand. The Medium TKN, built around the Army IFAK, consists of two main components – the outer MOLLE compatible pouch and the removable inner supply chassis. Inside are six rows of various sized heavy duty elastic which holds the Trauma Kit Refill or other medical supplies.

|  |                                 |
|--|---------------------------------|
| Weight   | 197g                            |
| Part #   | HW-TKN-MED-EMPTY                |
| Size   | 2.5" x 4.5" x 7"                |
| Features   |                                 |
| Patented Helium Whisper® MOLLE Attachment System | Single row of MOLLE on top flap |
| Materials  |                                 |
| 500D Solution Dyed CORDURA®                      | TEX 70 bonded nylon thread      |
| ULTRAcamp® laminate                              | Ten-Speed® Elastic              |



\*Other colors available upon request \*\*Accessories not included

MADE IN USA



877.430.2583  
Cage Code: 3X9S8  
DUNS Number: 152028291  
Sales@BlueForceGear.com  
166 Pine Barren Rd., Pooler, GA 31322



Specifications are approximate and subject to change without notice. ©2022 Blue Force Gear, Inc. All Rights Reserved.

**LIFETIME WARRANTY**

100% Berry Compliant